

pediatric seating

Special Tomato® seating system

Soft-Touch™ sitter seat

- latex-free, seamless material is impermeable to fluids and is peel/tear resistant
- anti-microbial protected
- mobile, stationary or wedge base sold separately
- soft contoured seat with anti-thrust seating surface
- built in abductor
- contoured head and lateral supports
- adjustable 5-point harness
- attachment straps to secure to standard chairs
- specify color: brown (BRN) or aqua (TEL)



sitter seat

30-3400 BRN/TEL	size 1	420.00
30-3401 BRN/TEL	size 2	552.50
30-3402 BRN/TEL	size 3	702.50
30-3403 BRN/TEL	size 4	777.50
30-3404 BRN/TEL	size 5	1,185.00

sitter with wedge (floor sitter)

- wedge for sitters sizes 1, 2, or 3 only
- straps secure the sitter to the wedge
- position the sitter on the wedge at various angles to provide the perfect support
- specify color of sitter: brown (BRN) or aqua (TEL)

wedge and sitter seat



sitter with wedge (floor sitter)

30-3460 BRN/TEL	sitter seat and wedge size 1	537.50
30-3461 BRN/TEL	sitter seat and wedge size 2	637.50
30-3462 BRN/TEL	sitter seat and wedge size 3	782.50

wedge only

30-3430 BRN	wedge only (brown only)	342.50
-------------	-------------------------	--------

sitter with stationary base

- includes sitter and stationary base
- base raises sitter off the floor
- position the sitter on the base 0-25° tilt
- specify color of sitter: brown (BRN) or aqua (TEL)



seat and stationary base

sitter with stationary base

30-3410 BRN/TEL	stationary sitter size 1	685.00
30-3411 BRN/TEL	stationary sitter size 2	822.50
30-3412 BRN/TEL	stationary sitter size 3	922.50
30-3413 BRN/TEL	stationary sitter size 4	1,000.00
30-3414 BRN/TEL	stationary sitter size 5	1,567.50

stationary base only

30-3415	stationary base for sizes 1-3	550.00
30-3416	stationary base for sizes 4-5	680.00

sitter with mobile base

- includes sitter and mobile base
- base has 4 swivel casters with foot operated rear locks (2)
- 25° tilt design
- specify color of sitter: brown (BRN) or aqua (TEL)



seat and mobile base

sitter with mobile base

30-3420 BRN/TEL	mobile sitter size 1	932.50
30-3421 BRN/TEL	mobile sitter size 2	1,040.00
30-3422 BRN/TEL	mobile sitter size 3	1,207.50
30-3423 BRN/TEL	mobile sitter size 4	1,332.50
30-3424 BRN/TEL	mobile sitter size 5	1,877.50

mobile base only

30-3425	mobile base for sizes 1-3	675.00
30-3426	mobile base for sizes 4-5	827.50

sizing chart	user weight	user height	inside depth	inside width	seat back height
Size 1	20-30 lb	30-35"	7"	8-10"	21"
Size 2	30-60 lb	35-48"	9½"	10-12"	26"
Size 3	40-80 lb	40-56"	11½"	12-14"	29"
Size 4	50-110 lb	50-65"	13½"	14-16"	35"
Size 5	90-200 lb	54-74"	16"	16-18"	39"