

car seats and strollers

Columbia® TheraPedic™ IPS™ car seat

- safely transport children with special needs
- shell side and positioning pads provide lateral support

standard features:

- high profile sides for optimum position and support
- vacuform ABS plastic shell and chromed steel side tubes
- durable, machine-washable, soft plush cover
- 5-point harness system
- crash and fire code 302 tested
- certified for use in aircraft
- available colors: blue (B) and rose (RS)
- child seat inside dimensions: 13"W at back x 10"D x 29"H; small adult seat inside dimensions: 17½"W at back x 12"D x 33"H

child seat



adult seat



31-1010	B/RS	IPS™ 2000 child seat: 20-102 lb, up 60"	1,095.00
31-1016	B/RS	IPS™ 2500 small adult seat: 40-130 lb, 54-66"	1,835.00

Columbia® TheraPedic™ IPS™ stroller transport

- lightweight folding strollers for use with TheraPedic™ IPS™ car seats
- ideal alternative to a traditional positioning wheelchair

standard features:

- five tilt angles: 20, 25, 30, 35, or 40 degrees
- height adjustable footrest
- large rear wheels for a smooth ride on any surface
- foldable frame for convenient storage and transport
- seat locks to safely secure car seat to mobility base
- open dimensions:
IPS™ 2100- 49"L x 22½"W x 36¼"H
IPS™ 2520- 54½"L x 25½"W x 37¾"H
IPS™ 2300- 36½"L x 25"W x 47½"H

small seat with stroller base



large seat with stroller base



umbrella stroller



stroller transport

31-1011	IPS™ 2100 stroller for small seat, 25 lb	3,195.00
31-1017	IPS™ 2520 stroller for large seat, 25 lb	3,195.00
31-1013	IPS™ 2300 umbrella stroller (seat optional), 22 lb	2,545.00
seat with stroller transport		
31-1012	small seat + IPS™ 2100 base	4,290.00
31-1018	large seat + IPS™ 2520 base	5,030.00
31-1014	small seat + IPS™ 2300 umbrella stroller	3,640.00

Trotter™ mobility chair

- lightweight and portable positioning transport chair
- adjustable seat angle, seat depth, back angle and swing-away footrest
- 5-point harness with reversible padded strap covers
- folds for easy transport and storage
- 8" front casters; 10" rear locking casters
- crash tested for transit
- black removable, washable fabric

chair		seat width	seat depth	seat-to-floor	back height	weight capacity	
31-1200	small	11 ¾"	11 ½" - 13 ½"	22"	20"	75 lb	1,100.00
31-1201	medium	13 ¾"	13" - 15"	22"	23"	100 lb	1,100.00
31-1202	large	15 ½"	13 ½" - 15 ½"	22"	24"	175 lb	1,150.00
31-1203	x-large	17 ½"	15" - 17"	22 ½"	25"	250 lb	1,150.00

accessories

31-1210	tray	220.00
31-1211	headrest extension	80.00
31-1212	padded head wings	140.00
31-1213	lateral supports	122.50
31-1214	foot/ankle positioner	102.50
31-1215	torso vest	217.50
31-1216	canopy	180.00
31-1217	bus transit tie-downs	140.00
31-1218	utility bag	65.00

