

Baseline[®] 5-position hydraulic pinch gauge



use without paddle to accommodate small hands (same pinch width as standard pinch gauges)

5-grip width can accommodate large hands

- easy adjustable paddle to accommodate to any size hand - other pinch gauges only have one pinch width
- user can't see results during testing
- 5 pinch positions permits additional tests for tip, palmar and key tests
 - 5-position pinch test
 - rapid exchange test
- pinch width without paddle is the same width as standard pinch gauge (mechanical and hydraulic)
- can now provide same test results as hydraulic hand dynamometers



12-0480	5-position clinical pinch gauge	250.00
12-0480P	additional paddle for 5-position clinical pinch gauge	25.00

Baseline[®] mechanical pinch gauge

- measure tip, key and palmar pinch strength in lb and kg
- accurate and repeatable measurements
- results are consistent with published Mathiowetz studies
- indicator remains at the maximum reading until reset
- includes protective case



for orthopedic and sports medicine

12-0200	30 lb (13.6 kg), blue	200.00
12-0201	60 lb (27 kg), red	200.00

for weak or damaged hand

12-0203	10 lb (4.5 kg), silver	200.00
12-0202	2 lb (0.9 kg), gold	200.00

Jamar[®] hydraulic pinch gauge



- 50 lb (23 kg) capacity
- accurate and repeatable pinch strength readings (tip, key and palmar)
- maximum strength indicator remains after each reading until reset
- dual scale shows lb and kg

12-0601	Jamar [®] hydraulic pinch	340.00
---------	------------------------------------	--------

Psytech[™] finger flexion/extension gauge



- lightweight plastic gauge measures finger flexion and extension
- separate scales record each measurement in both lb and kg
- maximum reading remains until reset
- flexion capacity is 20 lb; extension capacity is 5 lb

12-0230	flexion extension gauge	125.00
---------	-------------------------	--------

electronic pinch gauge



- 50 lb (23 kg) capacity
- easy-to-read LCD display
- includes button to zero last reading stored and two toggle buttons to accurately store up to 99 readings
- wrist strap and one 3V battery included
- 10"L x 5"W x 2"D (25 x 13 x 5cm)

12-0607	Jamar [®] electronic pinch gauge	252.50
12-0475	Baseline [®] electronic pinch gauge	200.00