


work hardening

lifting boxes



- test functional capacity of an individual's material handling abilities within work hardening programs
- with removable weight holders, hinged lid with stop mechanism

		L" x W" x H"	empty weight (lb)	weight cap (lb)		
A	55-1010	14 x 14 x 17	21½	100	standard	185.00
B	55-1011	8 x 8 x 12	9	25	small	135.00
C	55-1012	10½ x 10½ x 14½	14	50	medium	155.00
D	55-1013	13 x 13 x 16¼	19	80	large	165.00
E	55-1014	24 x 12 x 7¼	19	80	double	185.00
	55-1015	4-piece SET: 2 [55-1010], 1 [55-1011] & 1 [55-1012]				625.00

heavy duty sleds



- for testing, evaluating and conditioning individuals in work hardening and industrial rehab programs
- simulate everyday work motions
- control resistance and strengthening treatment for shoulders, back, hips, thighs and calves by adding weights
- anti-friction runners for a variety of surfaces
- sled specs: 36"L x 18"W; handle adjusts from 32" - 44"H; foam covered grips; 4 removable weight holders; 250 lb capacity
- remove weight holders and insert box to use with sled
- accessory box specs: 30"L x 15"W x 6"H; inside dimensions 14" x 22½"; removable rubber mat

F	55-1020	set: weight sled and accessory box	775.00
G	55-1021	weight sled only	635.00
H	55-1022	accessory box only	135.00

work device and push-pull weighted sled



- combination lifting box and push cart for work simulation, conditioning and evaluation
- strap, foam covered hand grips and hand slots for lifting
- cart promotes push/pull activities
- box is 13½"L x 13½"W x 12"H; cart is 17"L x 14"W x 4¾"H with 2" casters; T-handle adjusts 32-44" H; removable weight holder with weight keeper; 10 lb empty, 80 lb capacity
- companion sled can be used alone or with work device
- straight handle adjusts 36-48" H; interchangeable box; sled measures 16"L x 13½"W x 6"H; anti-friction runners

all accessories

I	55-1023	accessory box only	135.00
K	55-1030	cart with T-handle and accessory box	290.00
	55-1032	cart with straight handle and accessory box	310.00
J	55-1031	cart with T-handle	245.00
	55-1033	cart with straight handle	245.00
	55-1034	sled with T-handle and accessory box	475.00
	55-1036	sled with straight handle and accessory box	475.00
	55-1035	sled with T-handle	265.00
L	55-1037	sled with straight handle	250.00
K&L	55-1038	sled and cart, box, both handles	550.00

use Baseline® push-pull dynamometer to generate quantitative measurements

see pages 248-249 for more information and accessories

hydraulic analog (dial) readout

12-0392	50 lb / 23 kg	600.00
12-0393	100 lb / 45 kg	675.00
12-0394	250 lb / 113 kg	750.00
12-0388	500 lb / 226 kg	825.00

hydraulic digital (LCD) readout

12-0397	50 lb / 23 kg	900.00
12-0398	100 lb / 45 kg	975.00
12-0399	250 lb / 113 kg	1,050.00
12-0387	500 lb / 226 kg	1,125.00

electronic digital (LCD) readout

12-0340	50 lb / 23 kg	900.00
12-0341	100 lb / 45 kg	1,000.00
12-0342	250 lb / 113 kg	1,100.00
12-0343	500 lb / 226 kg	1,200.00

