

shoveling station



- multiple shelf heights allow for a varied range of arm, shoulder, leg and back motions which simulate repetitive on-the-job duties
- 3 angled shelves
- shovel included
- 72"L x 36"W x 60"H

55-1000 shovel station 765.00


work station



- wall mounted unit to improve manual dexterity
- height adjustable feature allows a variety of working postures
- widely spread pattern of bolt holes with a variety of fasteners
- 48"L x 10"W

55-1052 work station 305.00


hand assembly device



- table top device includes 102 parts for activities using bolts, washers and nuts that require either hand tightening or wrench
- hardware included
- 24"L x 8"W x 12"H

55-1053 assembly device 285.00

table top electrical assembly center



- table top center simulates job tasks involving assembly activities
- promotes upper extremities ROM as well as working in confined spaces
- basic electrical assembly supplies included (not for live current)

55-1054 assembly center 550.00

adjustable height shelf assembly and sorting bin

sorting bin

- encourages cognitive and perceptual skills
- used to develop ROM
- 4 adjustable panels on top; 6 fixed-sized bins on bottom
- 29½"L x 13"W x 16"H

adjustable height shelf and frame

- height adjustable work surface to evaluate and develop ROM, strength and proper lifting techniques
- 36"L x 36"W base; 34½"L x 12"W top shelf; 17¼"W x 29¾"L adjustable; shelf height from 6"-48"H in 6" increments



55-1040 shelf assembly w/sorting bin 750.00
55-1042 sorting bin only 165.00
55-1041 adjust. assembly only 585.00

deluxe electrical assembly



- stand alone center simulates job tasks and assembly activities
- encourages whole-body ROM and working in confined spaces
- basic electrical supplies included (not for live current)
- 22½"W x 10"D x 72"H; 28"W x 20"D base

55-1055 electrical assembly center 1,025.00

lifting/bolt center



- upper extremity ROM assembly box and work capacity box all in one
- for vocational rehab and working capacity evaluation
- includes removable weight holder, weight keeper and a variety of fasteners
- empty weight 11½ lb; 80 lb capacity
- 24"L x 13½"W x 12"H

55-1050 lifting/bolt center

230.00