


CanDo® Multi-Axial platform system

Multi-Axial Platform System offers adjustable ROM (range-of-motion) and weight training to maximize lower body rehabilitation. Includes reversible board (1 side for left foot and 1 side for right foot). 5-ball set (color-coded yellow, red, green, blue and black from easiest to most challenging). Balls screw into board to provide vestibular challenge. Comes standard with 2 weight rods, 10 lbs of weight plates and a wall rack. Optional motion limiter is used to limit ROM.



professional system can be hung on the wall for easy storage!

instability balls can be used for both Multi-Axial Platform System and MVP® systems



Professional System (pictured above)

- A two-sided professional Multi-Axial board
- B weight rack with five color-coded instability balls: yellow, red, green, blue and black
- C 2.5 lb weight plates (4 each)
- D motion limiter complete with clips and padding
- E weight rod (fits any of the five sockets on either side of board)

Multi-Axial Platform System - PRO

sets

10-1730	board with balls, rods	315.00
10-1731	board with balls, rods, weights, wall rack	425.00
10-1735	board with balls, rods, weights, storage tub	380.00

instability balls (hemispheres)

	challenge	height	each	
yellow	x-easy	1"	10-1760	10.00
red	easy	1.5"	10-1761	12.50
green	moderate	2"	10-1762	15.00
blue	hard	2.5"	10-1763	17.50
black	x-hard	3"	10-1764	20.00
yel red grn blu blk	5 ball set		10-1765	67.50
yel red grn blu blk	set with rack		10-1742	112.50

instability ball wall rack

10-1740	wall rack for 5 balls and board	50.00
---------	---------------------------------	-------

accessories

10-1732	Multi-Axial board only	275.00
10-1741	storage tub (not shown)	50.00
10-1733	motion limiter	50.00
10-0601	2.5 lb weight plate	5.00
10-0601-4	2.5 lb weight plate (4 each)	20.00
10-1770	weight rod (each)	7.50
50-1596B	Dycem® non-slip circular pad, 7½" dia, blue	15.00

Multi-Axial Platform System - Home



- for home use
- board is elliptical with off-centered balance pivot
- select board for left or right foot

adult model (over 120 lbs) - blue

10-1752	home board for right foot	45.00
10-1753	home board for left foot	45.00

child model (under 120 lbs) - yellow

10-1750	home board for right foot	40.00
10-1751	home board for left foot	40.00

non-skid accessory

50-1596B	Dycem® non-slip circular pad, 7½" dia, blue	15.00
----------	---	-------