


range-of-motion (ROM)

**increase upper extremity ROM, strength, coordination and endurance,
or improve cognitive and visual skills**

double curved shoulder arc



- allows customizing program for people with limited ROM
- moving the colored tabs from one side of the arc to the other achieves full upper extremity ROM and improved motor planning and visual tracking skills
- use while standing for balance
- 24½" high at peak, 7½" at valley
- includes 36"L x 3"W base and 24 colored tabs

10-1155 double curve 240.00

climbing board



- board and bar frame with 7 rubber-tipped extensions on each side spaced 4" apart
- 38" bar has handle grips
- can be attached to table by clamps or through base
- accessory kit contains pattern cards, plastic, numbered tabs and 100' nylon cord for additional activities
- 24½"L x 36"H

10-1152 board and bar 495.00
10-1153 accessory kit 135.00

horizontal ring tree



- for people with neurological disorders, burns and orthopedic problems
- includes 45 plastic rings (15 ea 1, 1½" and 2")
- 3 horizontal rods (12", 24" and 36"L)
- 36"L base

10-1158 horizontal tree 452.50

incline board



- molded-plastic board can be permanently mounted or clamped on a tabletop
- vary angle to provide different shoulder extensions
- bilateral or single handle
- handle piece: 2¾" x 5½" x 8½"
- weights can be added for increased resistance
- 24" base; 8½"W x 31"H

10-1161 incline board 427.50


vertical ring tree



- includes 6 plastic rods: (3 ea: 18"H and 13"H)
- rod diameters range from ¼" to 1"
- includes 12 rings for each rod
- 21¾"L x 14½"W x 1½"D

10-1156 vertical tree 275.00
10-1157 72 extra rings 87.50

semi-circular pegboard



- for people with hand injuries, neurological involvement or orthopedic injuries
- 3 tiers hold 43 pegs in total
- pegs come in 2 lengths and 3 diameters
- 26"L x 14"W x 2"H

10-1159 peg board 200.00

pipe-tree



- for bilateral coordination and problem solving
- includes 31 pieces of ½" PVC pipe
- box cover is used as base for construction
- box accommodates all the piping for convenient storage

10-0766 pipe tree 245.00