

paraffin heat therapy

Therabath®



- 6-9 lb paraffin capacity
- hospital grade paraffin warmer with lid, grill and guide
- inside unit dimensions: 11.4"L x 6.75"W x 5.8"D
- operating temperature of 130°F

11-1171	unit with 6 lb unscented wax	200.00
11-1170	unit with 6 lb wintergreen wax	200.00

Dickson® paraffin baths

- **Dickson PB-107:** 6 lb capacity, 12"L x 6"W x 6"H, stationary unit, timer controlled high-heat setting which returns to operating temperature after its cycle
- **Dickson PB-104:** 20 lb capacity, 14"L x 7"W x 9"H, includes removable mobile stand, timer-controlled sterilization and heated drain
- **Dickson PB-105:** 50 lb capacity, 25"L x 11.5"W x 10"H, includes removable mobile stand, timer-controlled sterilization and heated drain

A	11-1151	PB-107 with 6 lb paraffin	582.50
B	11-1161	mobile stand for PB-107	337.50
C	11-1152	PB-104 with stand with 20 lb paraffin	3,610.00
D	11-1153	PB-105 with stand with 50 lb paraffin	5,100.00
	11-1162	mobile stand for PB-104	325.00



Parabath®

- lightweight and portable for clinic or home use
- for the treatment of pain associated with arthritis, strains, sprains or joint stiffness
- also used to increase range-of-motion
- deep oversized tank (12"L x 6"W x 6"H) allows easy immersion of hand and wrist, foot and ankle, or elbow
- stainless steel tank with insulated outer housing to retain heat
- operating temperature is 125°-134°F

11-1110	unit with 6 lb unscented wax	212.50
11-1119	stand: fits 18"L x 10"W x 7"D unit	240.00



Whitehall® ParaTherapy®

- **Whitehall PT-06-s:** 6 lb capacity, hand paraffin bath 12"L x 6"W x 6"D
- **Whitehall PT-18-s:** 18 lb capacity, arm or hand paraffin bath 20"L x 11.5"W x 6"D

E	11-1148	PT-6-S with 6 lb paraffin	467.50
F	11-1149	PT-18-S with 18 lb paraffin	962.50

